Bioanalytics, Metabolomics and Pharmacokinetics Shared Resource (BMPK)

Director: Dr. James Mohler

Ceritinib in Heparinized Human Plasma

(Sensitivity: 1.00 ng/mL)

BMPK validated a highly sensitive liquid chromatographic tandem mass spectrometric assay (LC-MS/ MS) for the analysis of ceritinib in heparinized human plasma. Ceritinib (ZYKADIA[™]) has been approved by the FDA for the treatment of adult patients with anaplastic lymphoma kinase (ALK) positive metastatic non-small cell lung cancer (NSCLC) who have previously been treated with crizotinib. Compared to crizotinib, ceritinib is approximately 20-fold more potent and is more selective for ALK. Clinical studies of ceritinib co-administered with CYP3A inhibitors (e.g., ketoconazole) or inducers (e.g., rifampin) have resulted in changes in the overall exposure of ceritinib and should be avoided whenever possible. The validated method has been used to support a clinical trial conducted at Roswell Park Comprehensive Cancer Center entitled "Phase I Study of Ceritinib (LDK378), a Novel ALK Inhibitor, in Combination with Gemcitabine-Based Chemotherapy in Patients with Advanced Solid Tumors".

Specifications and Validation Performance			HN CH ₃
Matrix (Anticoagulant):	Human Plasma (Lithium Heparin)		
Required Volume:	100 μL		
Preparation Procedure:	Protein Precipitation		ĊH ₃
HPLC Column:	C18		Ceritinib
Mobile Phase:	Acetonitrile with Acetic Acid		Formula: C₂ଃH₃₀CIN₅O₃S MW: 558.1433 g/mol
Flow Rate:	300 μL/min		2000, 2000
Detection Type:	Tandem Mass Spectral (MS/MS)		1890
	4 00 500	<i></i>	1690 1800 1400 1400
Calibration Ranges:	1.00 - 500 ng/mL		1200 1200
Calibrator Accuracy: Calibrator Precision:	100% (97.8 - 102%; n=5) 1.79% CV (0.688 - 5.20%; n=5)		2 1008 1006 1 1008 800
Calibrator Precision:	1.79% CV (0.888 - 5.20%, 11-5)		600
QC Concentrations:	3.00, 15.0 and 375 ng/mL		400 400 200 200
QC Accuracy:	97.9% (95.6 - 101%; n=18)		
QC Precision:	4.91% CV (3.01 - 7.44%; n=18)		LLOQ S/N = 16.5 Matrix Blank
	Human Pharn	nacokinetic Parameters of Ceriti	nib ^{1,2}
Recommended Dosing		450 mg single oral daily administration without food; steady- state is achieved after ~15 days of daily dosing	
Maximum Tolerated Dose (MTD)		750 mg single oral daily administration	
Bioavailability		Absolute bioavailability unknown; systemic exposure increased when taken with food, which may increase adverse events	
Active Metabolites		None	
Metabolism		~82% of the drug circulates in plasma unchanged (~92% is elimi- nated in feces with ~68% as unchanged parent compound)	
Plasma Protein Binding		97.2% in humans; independent of drug concentration	
Maximum Plasma Concentration (C _{max})		Occurs 4-6 hours after dosing; increases dose proportionally over a 50 to 750 mg oral dose range	
Plasma Terminal Half-Life (t _{1/2})		41 hours with nonlinear PK over t	time
¹ Ceritinib (LDK378) Investigator's Bro		<u>^</u>	

BMPK offers a wide range of bioanalytical and PK/PD modeling services to assist investigators in their basic research, preclinical, and clinical study objectives. For information on services and pricing, contact Wenjuan Zha, Ph.D., Associate Director at (716) 845-3258 or Wenjuan.Zha@RoswellPark.org.

