HO

но

Bioanalytics, Metabolomics and Pharmacokinetics Shared Resource (BMPK)

Estrogens in Human Serum

(Sensitivity: 2.00 pg/mL for E1 and 1.00 pg/mL for E2)

BMPK has validated a highly sensitive liquid chromatographic tandem mass spectral assay (LC-MS/MS) for the analysis of estrone (E1) and 17β-estradiol (E2). These steroidal sex hormones are secreted in both males and females and exhibit a wide range of physiological activity from development and regulation of the female reproductive system and secondary sex characteristics to gene regulation (genomic effects) and cell signaling (epigenomic effects) via estrogen receptors ER α and ER β . The validated method was applied to the analysis of 823 serum samples obtained during a 5 year flaxseed study of postmenopausal women (performance data shown below). The assay was also adapted to the analysis of heparinized human plasma and breast tissue (benign and malignant) samples examining the mechanistic response of tamoxifen in a breast cancer study.

O CH ₃]]	Specifications and Performance	
Н	Matrix:	Human Serum
	Required Volume:	600 μL
но	Preparation Procedure:	Solid Phase Extraction
Estrone (E1)	HPLC Column:	C18⁺
	Mobile Phase:	Acetonitrile with Acetic Acid
	Flow Rate:	300 μL/min
H H	Detection Type:	Tandem Mass Spectral (MS/MS)
HHH	Calibration Range:	2.00 to 250 pg/mL for E1
но		1.00 to 250 pg/mL for E2
17β-Estradiol (E2)	QC Concentrations:	6.00, 30.0 and 180 pg/mL for E1
		3.00, 30.0 and 180 pg/mL for E2
	Calibrator Accuracy:	100% (97.3 - 106%; n=14) for E1
17β-E2		100% (97.4 - 102%; n=13) for E2
™ 17a-E2	Calibrator Precision:	3.26% CV (1.70 - 5.06%; n=14) for E1
1/0°E2		3.17% CV (1.74 - 5.36%; n=13) for E2
	QC Accuracy:	92.8% (91.1 - 94.6%; n=40) for E1
MA MA MA		96.9% (96.2 - 98.3%; n=37) for E2
NA INA INA INA INA INA INA INA INA INA I	QC Precision:	5.42% CV (4.06 - 7.67%; n=40) for E1
		4.78% CV (2.92 - 6.54%; n=37) for E2

BMPK offers a wide range of bioanalytical and PK/PD modeling services to assist investigators in their basic research, preclinical, and clinical study objectives. For information on services and pricing, contact Joshua Prey, MS, Research Project Administrator at (716) 845-3313 or Joshua.Prey@RoswellPark.org.

